

# BDL

Bridge Dental Laboratory Ltd

## ***SNOREGUARDS***



**BRIDGE DENTAL LABORATORY**

Info@bridgedentallab.com

01202 486113

www.bridgedentallab.co.uk

# MSI



Colour for photo purposes only, product will be clear.

Slim & low in profile, the MSi is designed with patient comfort and compliance in mind. It is also the only system with a dedicated protrusive registration (BiteRight) which eliminates the usual guesswork associated with wax bite provision.

- The MSi hook is extremely strong and all component parts can be replaced
- Extremely strong and hard wearing - withstands heavy bruxism
- Provides an impressive amount of protrusive adjustment - 11mm
- High level of comfort and patient compliance
- No plaque traps - all accessible parts
- Reliable first time fit due to the BiteRight jig and the large amount of adjustment
- Easy for patients to use and self-adjust - reduced need for extra appointments
- All parts are replaceable to ensure long life
- Allows normal breathing - with mouth open or closed.

If you have any questions about the MSI appliance, please get in touch  
01202 486113  
[info@bridgedentallab.com](mailto:info@bridgedentallab.com)

# SILENSOR



The Silensor-sl consists of one splint for the upper jaw and one splint for the lower jaw. The lower jaw is either held in a predetermined or advanced position by two connectors that are fixed laterally to the splint.

The Silensor-sl counteracts the narrowing of the respiratory tracts. As the velocity of the air decreases so do the noise-generating vibrations of the soft tissues. The Silensor-sl allows jaw movements so is comfortable to wear, but prevents the lower jaw from falling back so is an effective snoreguard.

Advantages include:

- improved patient comfort
- proven effectiveness in treating the symptoms of OSA
- completely metal-free
- easily adjustable

If you have any questions about the SILENSOR appliance,  
please get in touch

01202 486113

[info@bridgedentallab.com](mailto:info@bridgedentallab.com)

# FUSED



Our Fused snoreguard is a one-piece, non-adjustable appliance. The appliance is set with the mandible in a protruded position to maintain the airway and prevent it from dropping back during sleep to obstruct the airway.

Advantages include:

- More cost effective
- Easy to clean
- Great for partially edentulous patients
- Durable due to lack of components

If you have any questions about the FUSED appliance,  
please get in touch

01202 486113

[info@bridgedentallab.com](mailto:info@bridgedentallab.com)